

## 3T3 Neutral Red Uptake *in vitro* Phototoxicity Test

### Regulatory Acceptance

Phototoxicity (photo-irritation) is an acute light induced, non-immunologic skin response to a photo-reactive chemical. A thoroughly validated *in vitro* method is available for phototoxicity testing. The 3T3 NRU Phototoxicity test has been officially accepted by the EU Commission and the EU member states *Annex V to Directive 67/548/EEC on the Classification, Packaging and Labelling of Dangerous Substances: Testing Method B.41 Phototoxicity – in vitro 3T3 NRU Phototoxicity Test* and is also available as OECD guideline 432, *In Vitro 3T3 NRU Phototoxicity Test*, published November, 23, 2004.

### Test principle

The test is based on a comparison of the cytotoxicity of a chemical when tested in the presence and in the absence of exposure to a non-cytotoxic dose of UVA/vis light. The use of the permanent mouse fibroblast cell line Balb/c 3T3 is recommended in the standard protocol of the EU and the OECD guideline. Cytotoxicity is expressed as a concentration dependent reduction of the uptake of the vital dye, Neutral Red. In most cases data from this *in vitro* test may provide sufficient information for the preclinical assessment of the phototoxic potential of a drug product and thus *in vivo* non-clinical studies are normally not warranted.

Advantages of the 3T3 NRU Phototoxicity test include: regulatory approval for *in vitro* assessment of phototoxic potential, relatively simple technical procedures, excellent reproducibility of results and cost efficiency.

### Test procedure

Balb/c 3T3 cells are cultured in 96well plates to form monolayer. First, the 96well plates are incubated with different concentrations of the test substance. Some plates are exposed to a dose of light, while control plates are kept in the dark. Then, the assay medium is replaced with fresh culture medium and the cells are incubated for 18-20 hours. Afterwards, cell viability is determined with the vital dye neutral red. The optical density is determined using a microplate reader. Cell viability at the tested chemical concentration with and without light exposure is compared to determine the photo-irritation potential.

An UV lamp is used as a light source. A dose of UVA equivalent to 5 J/cm<sup>2</sup> is delivered to the 96well plate.

The IC<sub>50</sub> concentrations (i.e. the concentrations by which cellular viability is reduced by 50%) are determined for the plates incubated in the dark (Cytotoxicity) and under UVA (Phototoxicity), Then, the phototoxic potential of a test chemical or photo-irritancy factor (PIF) is determined by comparison of the IC<sub>50</sub> without and with UVA exposure:

$$PIF = \frac{IC_{50} (-UVA)}{IC_{50} (+UVA)}$$

In addition, the mean photo effect (MPE) is calculated. The MPE is based on the comparison of the complete concentration response curves. It is defined as the weighted average across a representative set of photo effect values. The MPE is a useful measure of phototoxic potential for cases in which to equally effective concentrations (IC<sub>50</sub>) in the dark (-UVA) and light (+UVA) cannot be determined.

Base on the guideline, the Phototoxicity potential of the test substance is interpreted as follows:

Mean Photo Effect (MPE)	Photo-irritancy Factor (PIF)	Phototoxic Potential
< 0.1	< 2	Non-phototoxic
> 0.1 and < 0.15	> 2 and < 5	Probably Phototoxic
> 0.15	> 5	Phototoxic