

**RabHCP-I Rabbit Hepatocytes Cryopreserved Plateable for Induction assays
Cell Specification – Certificate of Analysis (CoA)**

Lot RabH200330

Batch Release: June 29, 2020

Species: New Zealand white rabbit (*Oryctolagus cuniculus*)

Gender: male
Age: 2 months


Cryopreservation:

Date: March 30, 2020
Amount per vial: 7.0×10^6 cells


Thawing: n = 1

Post-thaw viability: 79.7 %
Post-thaw yield per vial: 3.8×10^6 cells
Recovery: 54 %

Phase contrast on day 1 after thawing
(24well plate)



Phase contrast on day 3 after thawing
(24well plate)



Viability test on orbital shaker (Eppendorf Thermomixer C, 1000 rpm at 37 °C with 0.5×10^6 cells in 0.5 ml HPM-Cryo):

Time (h)	0	1
Viability (%)	79.7	53.7

Recommended seeding density on collagen-coated plates:

24well plate – 300,000 cells/well // 96well plate – 50,000 cells/well

Culture in Human Hepatocyte Maintenance Medium (HHMM).

Use Corning Collagen coated plates.

CYP P450 activity in culture after thawing:

Ethoxyresorufin-O-deethylation:

Induction with 25 μ M β -Naphthoflavone

pmol/mL
24well: 90.7 ± 2.2

X-fold induction
2.3

Note: Yield, viability, recovery and activity assays were performed at PRIMACYT using PRIMACYT's manual for thawing, plating and culture of primary cryopreserved hepatocytes.

Store at -150 °C or in the vapour phase of LN₂

This product should be considered as potential biohazard. Only intended for *in vitro* use.

Issued by A. Ullrich

Verified by D. Runge